

NORTHLAND RED LIGHT THERAPY

OMNILUX CONTOUR LED FACE MASK

Frequently Asked Questions (FAQ)



What is the Omnilux Contour Mask?

The Omnilux Contour is a professional-grade LED light therapy mask designed to support healthier, brighter, firmer-looking skin using clinically researched red and near-infrared light wavelengths. It is completely non-invasive, painless, and requires no downtime.

How Does It Work?

The Omnilux Contour uses specific wavelengths of light that are absorbed by the skin and underlying tissues.

This process is called **photobiomodulation**, which means using light energy to stimulate normal cellular processes. The light helps support the skin's natural repair and rejuvenation processes.

Think of it like helping your skin cells function more efficiently.

Does It Use Heat?

No.

The treatment is powered by light, not heat.

While some people may notice slight warmth, the light itself is what creates the benefits. Omnilux LED therapy is considered a non-thermal treatment.

What Light Wavelengths Are Used?

Red Light (633nm)

Supports:

- ✓ Healthier-looking skin
- ✓ Improved skin tone and texture
- ✓ Cellular regeneration
- ✓ Increased circulation and tissue oxygenation
- ✓ Brighter and more radiant skin

Near-Infrared Light (830nm)

Supports:

- ✓ Collagen production
- ✓ Elastin production
- ✓ Deeper tissue processes
- ✓ Skin firmness and elasticity

The Contour mask combines both wavelengths in every treatment.

What Can the Omnilux Contour Help With?

Many clients use the Omnilux Contour to support:

- ✓ Overall skin health
 - ✓ Fine lines and wrinkles
 - ✓ Skin firmness
 - ✓ Skin elasticity
 - ✓ Uneven skin tone
 - ✓ Redness and inflammation
 - ✓ Dull or tired-looking skin
 - ✓ Recovery following certain cosmetic procedures
 - ✓ Skin brightness and hydration appearance
-

What Results Can I Expect?

Results are cumulative.

LED therapy is not a "one-and-done" treatment.

Most people notice improvements in stages:

Often Seen First

- ✓ Reduced redness
- ✓ Reduced inflammation
- ✓ Improved skin resilience

Then

- ✓ Improved texture
- ✓ Softer skin
- ✓ Brighter appearance

Later

- ✓ Firmer-looking skin
- ✓ Improved elasticity
- ✓ Softer appearance of fine lines and wrinkles

Consistent use provides the best results.

What Did Clinical Participants Report?

After using the Omnilux Contour three times per week for four weeks:

- ✓ 98% said their skin felt firmer
- ✓ 96% said fine lines appeared less visible
- ✓ 95% said their skin looked brighter and plumper

Participants also reported softer skin, improved tone, less redness, and reduced appearance of pigmentation.

How Long Is Each Session?

Treatment Time

 10 Minutes

The treatment is comfortable and hands-free.

Many clients choose to relax, meditate, or simply rest during treatment.

How Often Should I Use It?

For best results:

- ☒ 3-5 sessions per week
- ☒ 10 minutes per session
- ☒ For 4-6 weeks initially
- ☒ Then continue with maintenance sessions as desired

Consistency is more important than doing longer sessions.

Will It Hurt?

No.

The Omnilux Contour is painless, non-invasive, and requires no recovery time.

Many clients describe the experience as relaxing and comfortable.

Will I Feel Anything During Treatment?

Some people notice mild warmth.

Some people feel nothing at all.

Both experiences are completely normal.

Can I Wear Makeup During Treatment?

For best results, the mask should be used on clean, dry skin.

Removing makeup, sunscreen, and thick skincare products allows direct contact with the skin and helps optimize treatment. Based on provider guidance, skincare routines should always be considered when using LED therapy.

Can I Overuse the Mask?

More is not necessarily better.

Follow the recommended treatment schedule and focus on consistency.

Regular use over time produces the best results.

Is It Safe?

LED therapy is considered safe for most individuals.

However, you should consult with a healthcare provider before use if you:

- ⚠ Have a history of photosensitivity
- ⚠ Take medications that increase sensitivity to light
- ⚠ Use skincare products that may increase sensitivity to light
- ⚠ Have concerns about existing medical conditions

The provider training recommends discussing photosensitivity history, medications, and skincare products before beginning treatment.

Why Omnilux?

Omnilux is recognized as one of the most researched LED therapy brands available.

Highlights include:

- ✔ More than 20 years of continuous research
 - ✔ More than 40 peer-reviewed studies
 - ✔ Professional-grade technology
 - ✔ Trusted by dermatologists and aesthetic professionals worldwide
 - ✔ Flexible silicone design for a comfortable fit
-

Fun Fact

NASA research helped increase interest in LED light therapy when scientists observed that specific light wavelengths may support cellular repair and tissue health during space-related research.

Quick Talking Points

What is it? A professional LED light therapy mask designed to support healthier, brighter, firmer-looking skin.

How does it work? Specific wavelengths of red and near-infrared light help support normal cellular processes.

How long does it take? Just 10 minutes.

How often should I use it? 3-5 times per week for 4-6 weeks, then maintenance.

When will I see results? Results build over time. Most people notice reduced redness and improved skin texture before seeing improvements in firmness and fine lines.

Does it hurt? No. It's completely non-invasive, painless, and requires no downtime.

Northland Red Light Therapy

Restore • Recover • Recharge ✨

This treatment is intended to support skin wellness and appearance. It is not intended to diagnose, treat, cure, or prevent any disease. Kmkzuzf